

FOXN2

Reactivity:Human Mouse

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:47kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human FOXN2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HTLF;

Background:

This gene encodes a forkhead domain binding protein and may function in the transcriptional regulation of the human T-cell leukemia virus long terminal repeat.

To place an order, please [Click HERE](#).

Catalog #:A7185

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3344

Isotype:IgG

Swiss Prot:P32314

Purity:Affinity purification

For research use only.