

## GTF2H3

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:** 34kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human GTF2H3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

P34; BTF2; TFB4; TFIIH;

**Catalog #:** A7188

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 2967

**Isotype:** IgG

**Swiss Prot:** Q13889

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a member of the TFB4 family. The encoded protein is a subunit of the core-TFIIH basal transcription factor and localizes to the nucleus. The encoded protein is involved in RNA transcription by RNA polymerase II and nucleotide excision repair and associates with the Cdk-activating kinase complex. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 14.

*To place an order, please [Click HERE](#).*