

## KIAA1456

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:** 14kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human KIAA1456

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

TRM9L; C8orf79;

*To place an order, please [Click HERE](#).*

**Catalog #:** A7193

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 57604

**Isotype:** IgG

**Swiss Prot:** Q8N9K7

**Purity:** Affinity purification

For research use only.