

KLF9

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:27kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human KLF9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BTEB; BTEB1;

Background:

The protein encoded by this gene is a transcription factor that binds to GC box elements located in the promoter. Binding of the encoded protein to a single GC box inhibits mRNA expression while binding to tandemly repeated GC box elements activates transcription.

To place an order, please [Click HERE](#).

Catalog #:A7196

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:687

Isotype:IgG

Swiss Prot:Q13886

Purity:Affinity purification

For research use only.