

PECR

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:33kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PECR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TERP; DCRRP; PVIARL; HPDHASE; SDR29C1; HSA250303;

To place an order, please [Click HERE](#).

Catalog #:A7206

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:55825

Isotype:IgG

Swiss Prot:Q9BY49

Purity:Affinity purification

For research use only.