

SMCHD1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:226kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SMCHD1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A7214

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23347

Isotype:IgG

Swiss Prot:A6NHR9

Purity:Affinity purification

For research use only.

Background:

This gene encodes a protein which contains a hinge region domain found in members of the SMC (structural maintenance of chromosomes) protein family.

To place an order, please [Click HERE](#).