

TOB2

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW:37kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TOB2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TOB4; TOBL; TROB2;

Background:

TOB2 belongs to the TOB (see TOB1; MIM 605523)/BTG1 (MIM 109580) family of antiproliferative proteins, which are involved in the regulation of cell cycle progression.

To place an order, please [Click HERE](#).

Catalog #:A7223

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10766

Isotype:IgG

Swiss Prot:Q14106

Purity:Affinity purification

For research use only.