

TP53AIP1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 13kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TP53AIP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P53AIP1;

Catalog #: A7224

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 63970

Isotype: IgG

Swiss Prot: Q9HCN2

Purity: Affinity purification

For research use only.

Background:

This gene is specifically expressed in the thymus, and encodes a protein that is localized to the mitochondrion. The expression of this gene is inducible by p53, and it is thought to play an important role in mediating p53-dependent apoptosis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

To place an order, please [Click HERE](#).