

KDM5A

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 192kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human KDM5A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RBP2; RBBP2; RBBP-2;

Catalog #: A7238

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5927

Isotype: IgG

Swiss Prot: P29375

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the Jumonji, AT-rich interactive domain 1 (JARID1) histone demethylase protein family. The encoded protein plays a role in gene regulation through the histone code by specifically demethylating lysine 4 of histone H3. The encoded protein interacts with many other proteins, including retinoblastoma protein, and is implicated in the transcriptional regulation of Hox genes and cytokines. This gene may play a role in tumor progression.

To place an order, please [Click HERE](#).