

SETD6

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:51kDa/53kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human SETD6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A7242

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:79918

Isotype:IgG

Swiss Prot:Q8TBK2

Purity:Affinity purification

For research use only.