

BLVRB

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 22kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human BLVRB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FLR; BVRB; SDR43U1; HEL-S-10;

Catalog #: A7260

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 645

Isotype: IgG

Swiss Prot: P30043

Purity: Affinity purification

For research use only.

Background:

The final step in heme metabolism in mammals is catalyzed by the cytosolic biliverdin reductase enzymes A and B (EC 1.3.1.24).

To place an order, please [Click HERE](#).