

LMO4

Reactivity: Human Mouse Rat

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:50 - 1:100

Calculated MW: 18kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LMO4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A7268

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8543

Isotype: IgG

Swiss Prot: P61968

Purity: Affinity purification

For research use only.

Background:

This gene encodes a cysteine-rich protein that contains two LIM domains but lacks a DNA-binding homeodomain. The encoded protein may play a role as a transcriptional regulator or as an oncogene.

To place an order, please [Click HERE](#).