

DCD

Reactivity:Human Mouse

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:11kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DCD

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PIF; AIDD; DSEP; HCAP; DCD-1;

Catalog #:A7280

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:117159

Isotype:IgG

Swiss Prot:P81605

Purity:Affinity purification

For research use only.

Background:

This antimicrobial gene encodes a secreted protein that is subsequently processed into mature peptides of distinct biological activities. The C-terminal peptide is constitutively expressed in sweat and has antibacterial and antifungal activities. The N-terminal peptide, also known as diffusible survival evasion peptide, promotes neural cell survival under conditions of severe oxidative stress. A glycosylated form of the N-terminal peptide may be associated with cachexia (muscle wasting) in cancer patients. Alternative splicing results in multiple transcript variants encoding different isoforms.

To place an order, please [Click HERE](#).