

ING5

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 28kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ING5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

p28ING5;

Catalog #: A7288

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84289

Isotype: IgG

Swiss Prot: Q8WYH8

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is similar to ING1, a tumor suppressor protein that can interact with TP53, inhibit cell growth, and induce apoptosis. This protein contains a PHD-finger, which is a common motif in proteins involved in chromatin remodeling. This protein can bind TP53 and EP300/p300, a component of the histone acetyl transferase complex, suggesting its involvement in TP53-dependent regulatory pathway.

To place an order, please [Click HERE](#).