

## PDCD5

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:**14kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human PDCD5

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

TFAR19;

**Catalog #:**A7298

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**9141

**Isotype:**IgG

**Swiss Prot:**O14737

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a protein that is upregulated during apoptosis where it translocates rapidly from the cytoplasm to the nucleus. The encoded protein may be an important regulator of K(lysine) acetyltransferase 5 (a protein involved in transcription, DNA damage response and cell cycle control) by inhibiting its proteasome-dependent degradation. Pseudogenes have been identified on chromosomes 5 and 12

*To place an order, please [Click HERE](#).*