

## SETD5

---

**Reactivity:**Human Rat

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:**158kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human SETD5

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:**A7304

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**55209

**Isotype:**IgG

**Swiss Prot:**Q9C0A6

**Purity:**Affinity purification

For research use only.

**Background:**

This function of this gene has yet to be determined but mutations in this gene have been associated with autosomal dominant mental retardation-23. Alternative splicing results in multiple transcript variants encoding different isoforms.

**To place an order, please [Click HERE](#).**