

SETD8

Reactivity:Human Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:43kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SETD8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

SET8; KMT5A; SET07; PR-Set7;

To place an order, please [Click HERE](#).

Catalog #:A7305

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:387893

Isotype:IgG

Swiss Prot:Q9NQR1

Purity:Affinity purification

For research use only.