

## SGTA

---

**Reactivity:**Human

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:**34kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human SGTA

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

SGT; hSGT; alphaSGT;

**Catalog #:**A7306

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**6449

**Isotype:**IgG

**Swiss Prot:**O43765

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a protein which is capable of interacting with the major nonstructural protein of parvovirus H-1 and 70-kDa heat shock cognate protein; however, its function is not known. Since this transcript is expressed ubiquitously in various tissues, this protein may serve a housekeeping function.

*To place an order, please [Click HERE](#).*