

SMYD3

Reactivity:Human Mouse

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:49kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SMYD3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

KMT3E; ZMYND1; ZNFN3A1; bA74P14.1;

Catalog #:A7309

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:64754

Isotype:IgG

Swiss Prot:Q9H7B4

Purity:Affinity purification

For research use only.

Background:

This gene encodes a histone methyltransferase which functions in RNA polymerase II complexes by an interaction with a specific RNA helicase. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).