

WBSCR22

Reactivity: Human Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 32kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human WBSCR22

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

WBMT; MERM1; PP3381; HUSSY-3; HASJ4442;

Catalog #: A7317

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 114049

Isotype: IgG

Swiss Prot: O43709

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein containing a nuclear localization signal and an S-adenosyl-L-methionine binding motif typical of methyltransferases, suggesting that the encoded protein may act on DNA methylation. This gene is deleted in Williams syndrome, a multisystem developmental disorder caused by the deletion of contiguous genes at 7q11.23. Alternatively spliced transcript variants have been found.

To place an order, please [Click HERE](#).