

PAK2

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:58kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PAK2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PAK65; PAKgamma;

Catalog #:A7333

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5062

Isotype:IgG

Swiss Prot:Q13177

Purity:Affinity purification

For research use only.

Background:

The p21 activated kinases (PAK) are critical effectors that link Rho GTPases to cytoskeleton reorganization and nuclear signaling. The PAK proteins are a family of serine/threonine kinases that serve as targets for the small GTP binding proteins, CDC42 and RAC1, and have been implicated in a wide range of biological activities. The protein encoded by this gene is activated by proteolytic cleavage during caspase-mediated apoptosis, and may play a role in regulating the apoptotic events in the dying cell.

To place an order, please [Click HERE](#).