

## RAB5C

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 23kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human RAB5C

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

RABL; L1880; RAB5L; RAB5CL;

**Catalog #:** A7342

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 5878

**Isotype:** IgG

**Swiss Prot:** P51148

**Purity:** Affinity purification

For research use only.

**Background:**

Members of the Rab protein family are small GTPases of the Ras superfamily that are thought to ensure fidelity in the process of docking and/or fusion of vesicles with their correct acceptor compartment (Han et al., 1996)

*To place an order, please [Click HERE](#).*