

VLDLR

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:96kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human VLDLR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CARMQ1; CHRMQ1; VLDLRCH;

Catalog #:A7345

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7436

Isotype:IgG

Swiss Prot:P98155

Purity:Affinity purification

For research use only.

Background:

The low density lipoprotein receptor (LDLR) gene family consists of cell surface proteins involved in receptor-mediated endocytosis of specific ligands. This gene encodes a lipoprotein receptor that is a member of the LDLR family and plays important roles in VLDL-triglyceride metabolism and the reelin signaling pathway. Mutations in this gene cause VLDLR-associated cerebellar hypoplasia. Alternative splicing generates multiple transcript variants encoding distinct isoforms for this gene.

To place an order, please [Click HERE](#).