

MAP4K3

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:101kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human MAP4K3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GLK; MEKKK3; MEKKK 3; MAPKKKK3; RAB8IPL1;

Background:

This gene encodes a member of the mitogen-activated protein kinase kinase kinase family. The encoded protein activates key effectors in cell signalling, among them c-Jun. Alternatively spliced transcripts encoding multiple isoforms have been observed for this gene.

To place an order, please [Click HERE](#).

Catalog #:A7351

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8491

Isotype:IgG

Swiss Prot:Q81VH8

Purity:Affinity purification

For research use only.