

ZNF264

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 71kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ZNF264

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A7356

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9422

Isotype: IgG

Swiss Prot: O43296

Purity: Affinity purification

For research use only.

Background:

This gene encodes a zinc finger protein and belongs to the krueppel C2H2-type zinc-finger protein family. Zinc finger proteins are often localized in the nucleus, bind nucleic acids, and regulate transcription.

To place an order, please [Click HERE](#).