

ERGIC2

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 43kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ERGIC2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PTX1; CDA14; Erv41; cd002;

Catalog #: A7369

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51290

Isotype: IgG

Swiss Prot: Q96RQ1

Purity: Affinity purification

For research use only.

Background:

ERGIC2, or PTX1, is a ubiquitously expressed nuclear protein that is downregulated in prostate carcinoma (Kwok et al., 2001)

To place an order, please [Click HERE](#).