

## CAMK1G

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 53kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human CAMK1G

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

VWS1; CLICK3; CLICKIII; dJ272L16.1;

**Catalog #:** A7379

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 57172

**Isotype:** IgG

**Swiss Prot:** Q96NX5

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a protein similar to calcium/calmodulin dependent protein kinase, however, its exact function is not known.

*To place an order, please [Click HERE](#).*