

NLRC4

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:116kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human NLRC4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CLAN; IPAF; AIFEC; CLAN1; CLANA; CLANB; CLANC; CLAND; FCAS4; CARD12; CLR2.1;

Background:

This gene encodes a member of the caspase recruitment domain-containing NLR family. Family members play essential roles in innate immune response to a wide range of pathogenic organisms, tissue damage and other cellular stresses. Mutations in this gene result in autoinflammation with infantile enterocolitis. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).

Catalog #:A7382

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:58484

Isotype:IgG

Swiss Prot:Q9NPP4

Purity:Affinity purification

For research use only.