

NAA50

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:19kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human NAA50

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SAN; MAK3; NAT5; hSAN; NAT13; NAT5P; hNAT5; NAT13P;

To place an order, please [Click HERE](#).

Catalog #:A7387

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:80218

Isotype:IgG

Swiss Prot:Q9GZZ1

Purity:Affinity purification

For research use only.