

CSRNP2

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:60kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CSRNP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C12orf2; PPP1R72; TAIP-12; C12orf22; FAM130A1;

Background:

The protein encoded by this gene belongs to the CSRNP family of nuclear proteins that share conserved regions, including cysteine- and serine- rich regions, a basic domain, a transcriptional activation domain, and bind the sequence 'AGAGTG', thus have the hallmark of transcription factors. Studies in mice suggest that these genes may have redundant functions. Alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).

Catalog #:A7388

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:81566

Isotype:IgG

Swiss Prot:Q9H175

Purity:Affinity purification

For research use only.