

SETDB2

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:82kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SETDB2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CLLD8; CLLL8; KMT1F; C13orf4;

Background:

Proteins that contain a SET domain, such as SETDB2, modulate gene expression epigenetically through histone H3 (see MIM 601128) methylation. SETDB2 is likely a histone H3 methyltransferase, as it contains both the active site and flanking cysteine residues required for catalytic activity (Zhang et al., 2003)

To place an order, please [Click HERE](#).

Catalog #:A7391

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:83852

Isotype:IgG

Swiss Prot:Q96T68

Purity:Affinity purification

For research use only.