

RIPK1

Reactivity:Human

Tested applications:WB IP

Recommended Dilution:WB 1:500 - 1:2000 IP 1:20 - 1:50

Calculated MW:76kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RIPK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RIP; RIP1;

To place an order, please [Click HERE](#).

Catalog #:A7414

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8737

Isotype:IgG

Swiss Prot:Q8NFZ0

Purity:Affinity purification

For research use only.