

ACAA1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:44kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ACAA1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ACAA; THIO; PTHIO;

Background:

This gene encodes an enzyme operative in the beta-oxidation system of the peroxisomes. Deficiency of this enzyme leads to pseudo-Zellweger syndrome. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).

Catalog #:A7422

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:30

Isotype:IgG

Swiss Prot:P09110

Purity:Affinity purification

For research use only.