

## ANXA11

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

**Calculated MW:**54kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human ANXA11

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ANX11; CAP50;

**Catalog #:**A7423

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**311

**Isotype:**IgG

**Swiss Prot:**P50995

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a member of the annexin family, a group of calcium-dependent phospholipid-binding proteins. Annexins have unique N-terminal domains and conserved C-terminal domains, which contain the calcium-dependent phospholipid-binding sites. The encoded protein is a 56-kD antigen recognized by sera from patients with various autoimmune diseases. Several transcript variants encoding two different isoforms have been identified.

*To place an order, please [Click HERE](#).*