

BLK

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:58kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human BLK

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MODY11;

Background:

This gene encodes a nonreceptor tyrosine-kinase of the src family of proto-oncogenes that are typically involved in cell proliferation and differentiation. The protein has a role in B-cell receptor signaling and B-cell development. The protein also stimulates insulin synthesis and secretion in response to glucose and enhances the expression of several pancreatic beta-cell transcription factors.

To place an order, please [Click HERE](#).

Catalog #:A7427

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:640

Isotype:IgG

Swiss Prot:P51451

Purity:Affinity purification

For research use only.