

AKR1C4

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:37kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human AKR1C4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C11; CDR; DD4; CHDR; DD-4; HAKRA; 3-alpha-HSD;

Background:

This gene encodes a member of the aldo/keto reductase superfamily, which consists of more than 40 known enzymes and proteins. These enzymes catalyze the conversion of aldehydes and ketones to their corresponding alcohols by utilizing NADH and/or NADPH as cofactors. The enzymes display overlapping but distinct substrate specificity. This enzyme catalyzes the bioreduction of chlordecone, a toxic organochlorine pesticide, to chlordecone alcohol in liver. This gene shares high sequence identity with three other gene members and is clustered with those three genes at chromosome 10p15-p14.

To place an order, please [Click HERE](#).

Catalog #:A7430

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1109

Isotype:IgG

Swiss Prot:P17516

Purity:Affinity purification

For research use only.