

SPIB

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:29kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SPIB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SPI-B;

Catalog #:A7451

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6689

Isotype:IgG

Swiss Prot:Q01892

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a transcriptional activator that binds to the PU-box (5'-GAGGAA-3') and acts as a lymphoid-specific enhancer. Four transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).