

FGFBP1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:26kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human FGFBP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FGFBP; HBP17; FGF-BP; FGF-BP1; FGFBP-1;

Background:

This gene encodes a secreted fibroblast growth factor carrier protein. The encoded protein plays a critical role in cell proliferation, differentiation and migration by binding to fibroblast growth factors and potentiating their biological effects on target cells. The encoded protein may also play a role in tumor growth as an angiogenic switch molecule, and expression of this gene has been associated with several types of cancer including pancreatic and colorectal adenocarcinoma. A pseudogene of this gene is also located on the short arm of chromosome 4.

To place an order, please [Click HERE](#).

Catalog #:A7458

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9982

Isotype:IgG

Swiss Prot:Q14512

Purity:Affinity purification

For research use only.