

PCGF3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 28kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PCGF3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RNF3; DONG1; RNF3A;

Catalog #: A7459

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10336

Isotype: IgG

Swiss Prot: Q3KNV8

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene contains a C3HC4 type RING finger, which is a motif known to be involved in protein-protein interactions. The specific function of this protein has not yet been determined.

To place an order, please [Click HERE](#).