

PCGF3

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:28kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PCGF3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RNF3; DONG1; RNF3A;

Background:

The protein encoded by this gene contains a C3HC4 type RING finger, which is a motif known to be involved in protein-protein interactions. The specific function of this protein has not yet been determined.

To place an order, please [Click HERE](#).

Catalog #:A7459

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10336

Isotype:IgG

Swiss Prot:Q3KNV8

Purity:Affinity purification

For research use only.