

VENTX

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 28kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human VENTX

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NA88A; HPX42B; VENTX2;

Catalog #: A7467

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27287

Isotype: IgG

Swiss Prot: O95231

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the Vent family of homeodomain proteins. The encoded protein may function as a transcriptional repressor and be involved in mesodermal patterning and hemopoietic stem cell maintenance. Multiple pseudogenes exist for this gene. A transcribed pseudogene located on chromosome X may lead to antigen production in certain melanomas.

To place an order, please [Click HERE](#).