

RPGRIP1

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:147kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RPGRIP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LCA6; RGI1; RGRIP; CORD13; RPGRIP; RPGRIP1d;

Background:

This gene encodes a photoreceptor protein that interacts with retinitis pigmentosa GTPase regulator protein and is a key component of cone and rod photoreceptor cells. Mutations in this gene lead to autosomal recessive congenital blindness.

To place an order, please [Click HERE](#).

Catalog #:A7477

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:57096

Isotype:IgG

Swiss Prot:Q96KN7

Purity:Affinity purification

For research use only.