

ELF2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 64kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ELF2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

EU32; NERF; NERF-2; NERF-1A; NERF-1B; NERF-1a,b;

To place an order, please [Click HERE](#).

Catalog #: A7487

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1998

Isotype: IgG

Swiss Prot: Q15723

Purity: Affinity purification

For research use only.