

HCCS

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:31kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human HCCS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MLS; CCHL; MCOPS7;

Catalog #:A7490

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3052

Isotype:IgG

Swiss Prot:P53701

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is an enzyme that covalently links a heme group to the apoprotein of cytochrome c. Defects in this gene are a cause of microphthalmia syndromic type 7 (MCOPS7). Three transcript variants encoding the same protein have been found for this gene.

To place an order, please [Click HERE](#).