

## RAC3

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

**Calculated MW:**21kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human RAC3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

q

**Catalog #:**A7498

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**5881

**Isotype:**IgG

**Swiss Prot:**P60763

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is a GTPase which belongs to the RAS superfamily of small GTP-binding proteins. Members of this superfamily appear to regulate a diverse array of cellular events, including the control of cell growth, cytoskeletal reorganization, and the activation of protein kinases.

*To place an order, please [Click HERE](#).*