

RAB31

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 22kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RAB31

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Rab22B;

Catalog #: A7506

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11031

Isotype: IgG

Swiss Prot: Q13636

Purity: Affinity purification

For research use only.

Background:

Small GTP-binding proteins of the RAB family, such as RAB31, play essential roles in vesicle and granule targeting (Bao et al., 2002)

To place an order, please [Click HERE](#).