

TRIM29

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:66kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TRIM29

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ATDC;

Catalog #:A7508

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23650

Isotype:IgG

Swiss Prot:Q14134

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the TRIM protein family. It has multiple zinc finger motifs and a leucine zipper motif. It has been proposed to form homo- or heterodimers which are involved in nucleic acid binding. Thus, it may act as a transcriptional regulatory factor involved in carcinogenesis and/or differentiation. It may also function in the suppression of radiosensitivity since it is associated with ataxia telangiectasia phenotype.

To place an order, please [Click HERE](#).