

UBA6

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 118kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human UBA6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

E1-L2; MOP-4; UBE1L2;

Catalog #: A7511

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55236

Isotype: IgG

Swiss Prot: A0AVT1

Purity: Affinity purification

For research use only.

Background:

Modification of proteins with ubiquitin (UBB; MIM 191339) or ubiquitin-like proteins controls many signaling networks and requires a ubiquitin-activating enzyme (E1), a ubiquitin conjugating enzyme (E2), and a ubiquitin protein ligase (E3). UBE1L2 is an E1 enzyme that initiates the activation and conjugation of ubiquitin-like proteins (Jin et al., 2007)

To place an order, please [Click HERE](#).