

SESN2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 54kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SESN2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HI95; SES2; SEST2;

Catalog #: A7515

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 83667

Isotype: IgG

Swiss Prot: P58004

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the sestrin family of PA26-related proteins. The encoded protein may function in the regulation of cell growth and survival. This protein may be involved in cellular response to different stress conditions.

To place an order, please [Click HERE](#).