

MCC

Reactivity:Human

Tested applications:IHC

Recommended Dilution:IHC 1:50 - 1:200

Calculated MW:93kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human MCC

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MCC1;

Catalog #:A7529

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4163

Isotype:IgG

Swiss Prot:P23508

Purity:Affinity purification

For research use only.

Background:

This gene is a candidate colorectal tumor suppressor gene that is thought to negatively regulate cell cycle progression. The orthologous gene in the mouse expresses a phosphoprotein associated with the plasma membrane and membrane organelles, and overexpression of the mouse protein inhibits entry into S phase. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).